
**DENTSCARE LTDA** - Av. Edgar Nelson Meister, 474 - Distrito Industrial - 89219-501 - Joinville - SC - Brasil / Brazil. CNPJ / Tax ID: 05.106.945/0001-06 INDUSTRIA BRASILEIRA / BRAZILIAN INDUSTRY - Responsável Técnico: Friedrich Georg Mittelstädt CTRC: 13100147 - SC - ANVISA nº 80172319014 - Marca/ Brand: FGM® / Fabricado no Brasil / Made in Brazil

**NUNO FLORES** - Al. Bonifácio Lázaro Lozano, 3 – Piso 0 – C 2780-125 Oeiras / Portugal - Tel: (351) 21 4439292 - ec.rep@fgm.ind.br

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### Sistema Adesivo Auto-condicionante para Esmalte e Dentina

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that interact among each other amplifying the curing capacity of light emitted from light-curing units. Added to different materials, the system provides different advantages.  
 When incorporated to Ambar Universal APS's composition, APS increases the conversion degree in the hybrid layer, which enhances the adhesive's bond strength and increases mechanical properties of the adhesive film (greater cohesive strength). Now the product is colorless, which provides a remarkable visual advantage when compared to other adhesive systems. This avoids any type of interference when performing a restoration/cementation on anterior teeth.  
**Product Presentation**  
 Package containing a bottle with 5ml of the product.  
**Basic Composition**  
 Active Ingredients: MDP (10-methacryloyloxydecyl dihydrogen phosphate), methacrylate monomers, photoinitiators complex (APS), co-initiators, stabilizers.  
 Inactive Ingredients: inert load (silica particles) and vehicle (ethanol).  
**Product Indication**  
 - All classes of direct restoration with composites (class I, II, III, IV and V)  
 - Adhesive luting (in addition to resin cement) of prosthetic pieces (post/cores, crowns, onlays/inlays, veneers, etc.) made of fiberglass, ceromer, ceramic, resin and metal;  
 - Adhesive repairs in ceramic and composites;  
 - Use as metallic or ceramic primer.  
**Precautions and Contraindications**  
 1. Before using the bonding agent, make sure that the dentin surface is slightly shiny (moistened), but not wet;  
 2. The air used in the whole procedure must be oil and water free;  
 3. The recommended light-curing time for this product depends on the power of the light emitted by the equipment used. The recommended power intensity of the light-curing unit is above 450mW/cm² which guarantees the quality of the polymer generated and the quality of the bonding. Light-curing units must have its light intensity checked by radiometer frequently. It is recommended that the patient and the surgeon use protective goggles during the emission of blue light by the light-curing unit;  
 4. Do not apply excess of adhesive to avoid compromising the quality of adhesion and prevent any aesthetic and functional problems. Excessive application of air blasts should also be avoided to prevent air bubbles in the adhesive;  
 5. Do not leave the adhesive at the reach of patients and children;  
 6. The adhesive contains methacrylate monomers. Avoid using the product in patients with allergy to acrylates/ methacrylates. To reduce risks, minimize the exposure time of the dental tissues to the material when it is not polymerized;  
 7. Wear protective gloves when handling the product. If the product comes into accidental contact with eyes or oral tissues, wash the area abundantly with water;  
 8. Make sure to close the bottle lid after use;  
 9. In deep cavities protect the dentin-pulp complex with a base or liner material. Do not use material that contains Eugenol as a basis to avoid insufficient curing on the surface. The application of high concentration of H<sub>2</sub>O<sub>2</sub> should also be avoided because it

can interfere on the curing process of the adhesive;  
 10. There is no need for shaking the bottle before using the adhesive.  
**Side Effects**  
 The adhesive can cause sensitization reactions in patients who are sensitive to any components of the formula. In this case, the material should not be used.  
**Instructions for Use**  
 Before starting the treatment read carefully the Warnings, Precautions, Contraindications and possible Side Effects.  
 Rewetting the dentin for the application of the adhesive system is necessary. Use cotton, pellets moistened with distilled water or chlorhexidine 2%. By rewetting the dentin the enamel is also rewetted, which can be treated by using the moist technique.  
 a. Cementation of prefabricated post: Apply phosphoric acid gel 37% into the root canal and remaining tooth structure for 15 seconds. Wash the area carefully and remove the excess water with absorbent paper cones.  
**3. ADHESIVE APPLICATION:**  
 3.1 Direct and Indirect Restorations: Dispense Ambar Universal APS in a Dappen pot or directly onto a disposable micro applicator. Apply two layers of adhesive – one drop for each– on the slightly moistened tooth surface. The first layer should be applied vigorously by rubbing the adhesive microapplicator saturated with the product for 10 seconds. Next, the second layer of adhesive – with a new drop – is applied for another 10 seconds and then the area is air blasted gently for 10 seconds to evaporate the solvent. Light cure Ambar Universal APS with blue light for 10 seconds.  
 3.2 Cementation of prefabricated posts: With the aid of a disposable micro applicator, apply two layers – one drop for each of the product into the root canal vigorously rubbing the adhesive on the walls of the conduit for 10 seconds. Evaporate the solvent with gentle blasts of oil-free air for 10 seconds, and then light cure it for 20 seconds.  
 OBS: During the adhesive application for post cementation, avoid accumulation of the product in the apical region. Excess of adhesive in this region compromises the evaporation of the solvent and polymerization.  
 3.3 As metallic or non-etchable ceramics primer: After sandblasting the internal surface of the prosthetic pieces, wash and air dry it. With the aid of a disposable micro applicator, apply Ambar Universal APS under friction on the internal surface of the previously treated piece, leave it undisturbed for 15 seconds and then apply gentle air for 10 seconds. Do not light cure, then apply the resin cement on the internal surface of the piece.  
 3.4 As etchable ceramics primer: The application method of Ambar Universal APS is equivalent to that indicated on item 3.3, however instead of sandblasting the surface, it must be acid etched with fluoridric acid and receive silane agent before Ambar Universal APS application. For further instructions on how to apply the acid and silane please refer to the respective instruction manuals.  
 - The recommended cure time for this product depends on the power of the light emitted by the equipment used. The recommended power density of the light-curing unit is above 450mW/cm² which guarantees the quality of the polymer generated and also the quality of the bonding.  
**4. BONDING OF RESTORATIVE MATERIALS:**  
 Perform restoration or cementation according to specific instructions of the product selected for the procedure

37% only on enamel for 15 seconds. Wash the surface with water and dry the cavity until it is moist, not dehydrated.  
 b. Luting of prefabricated posts: there is no need for prior etching of dentin.  
 2.3 Total-etch method:  
 a. Direct and Indirect Restorations: Apply phosphoric acid 37% for 15 seconds on both enamel and dentin, according to the type of preparation. Wash the surface with water and dry the cavity so that the dentin does not get dehydrated but also there is no accumulation of water on the surface. If rewetting the dentin for the application of the adhesive system is necessary, use cotton, pellets moistened with distilled water or chlorhexidine 2%. By rewetting the dentin the enamel is also rewetted, which can be treated by using the moist technique.  
 a. Cementation of prefabricated post: Apply phosphoric acid gel 37% into the root canal and remaining tooth structure for 15 seconds. Wash the area carefully and remove the excess water with absorbent paper cones.  
**3. ADHESIVE APPLICATION:**  
 3.1 Direct and Indirect Restorations: Dispense Ambar Universal APS in a Dappen pot or directly onto a disposable micro applicator. Apply two layers of adhesive – one drop for each– on the slightly moistened tooth surface. The first layer should be applied vigorously by rubbing the adhesive microapplicator saturated with the product for 10 seconds. Next, the second layer of adhesive – with a new drop – is applied for another 10 seconds and then the area is air blasted gently for 10 seconds to evaporate the solvent. Light cure Ambar Universal APS with blue light for 10 seconds.  
 3.2 Cementation of prefabricated posts: With the aid of a disposable micro applicator, apply two layers – one drop for each of the product into the root canal vigorously rubbing the adhesive on the walls of the conduit for 10 seconds. Evaporate the solvent with gentle blasts of oil-free air for 10 seconds, and then light cure it for 20 seconds.  
 OBS: During the adhesive application for post cementation, avoid accumulation of the product in the apical region. Excess of adhesive in this region compromises the evaporation of the solvent and polymerization.  
 3.3 As metallic or non-etchable ceramics primer: After sandblasting the internal surface of the prosthetic pieces, wash and air dry it. With the aid of a disposable micro applicator, apply Ambar Universal APS under friction on the internal surface of the previously treated piece, leave it undisturbed for 15 seconds and then apply gentle air for 10 seconds. Do not light cure, then apply the resin cement on the internal surface of the piece.  
 3.4 As etchable ceramics primer: The application method of Ambar Universal APS is equivalent to that indicated on item 3.3, however instead of sandblasting the surface, it must be acid etched with fluoridric acid and receive silane agent before Ambar Universal APS application. For further instructions on how to apply the acid and silane please refer to the respective instruction manuals.  
 - The recommended cure time for this product depends on the power of the light emitted by the equipment used. The recommended power density of the light-curing unit is above 450mW/cm² which guarantees the quality of the polymer generated and also the quality of the bonding.  
**4. BONDING OF RESTORATIVE MATERIALS:**  
 Perform restoration or cementation according to specific instructions of the product selected for the procedure

## EN Instructions for use

### Self-Etch Adhesive System for Enamel and Dentin

**Professional Use Only**  
*Read all this information carefully before using the product. Keep it for later reference until total consumption of the product and/or until there is no more interaction between the product and its last patient.*

**Product description**  
 Ambar Universal APS is a single bottle self-etch light-curing adhesive system (Seventh Generation) responsible for the union between tooth structure (enamel and dentin) and restorative materials. The primer and bonding agent are combined in a single bottle. Prior use of acid etchant is optional on both enamel and dentin, it can be used, therefore, in three protocols: self-etch, selective-etch and total-etch.

It can be used for direct restorations made of composite resins and light-curing glass ionomer or for indirect procedures when associated with resin cement. The product contains MDP, which allows it to be used as a metallic and ceramic primer.  
 APS is the acronym for Advanced Polymerization System, and it consists of a combination of different photoinitiators

(composé resins and/or cements).

**Conservation and Storage**

- Ambar Universal APS must be stored in temperatures between 5°C to 25°C/ 41°F – 77°F.

· Protect product from external contamination;

· The bottle containing the product must be closed carefully after use in order to prevent solvent evaporation.

· Do not expose the product to intense light or high temperatures.

· Do not use the product after expiration date.

· For cleaning, Ambar Universal APS can be removed with alcohol.

- Ambar Universal APS expires after 2 years from date of manufacture which is described on the box of the product.

**Warnings:**

Do not use the product after expiration date. In order to dispose the product properly follow your country’s legislation. Do not reuse empty packages. Keep out of the reach of children.

*This material was manufactured for dental use only and should be handled in accordance with the instructions. The manufacturer is not liable for damage caused by other uses or by mishandling. In addition, the user is obliged to prove, if this material is compatible with the desired use before using it and under his responsibility, especially when such use is not indicated in this instruction. Overview of data does not constitute any kind of guarantee and, therefore, has no binding.*

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