PRODUCTS 2019





Precision impressions even in extreme situations.

Elastomeric light body, medium body, heavy body and kneadable A-silicone impression materials for taking one-step and two-step impressions with modern curing characteristics.

Product properties

2

- Conventional A-silicones with maximum elastic recovery
- Most current setting characteristics
- "Genuine" putty from the 5:1 cartridge system or manually kneadable
- Diverse light-bodies with pronounced flowability and optimum hydrophilicity

User advantages

Economic material use, precisely dispensable

Over 30 years of

innovation quality and experience

- Reduced risk of errors due to quick-setting characteristics
- ✓ Highly precise final impression phase
- Different viscosities for many impression techniques and indications

Benefits for the dentist

- Maximum reliability and very easy handling when taking and removing impressions
- Elastomeric and stable storage, thixotropic, hydrophilic and flowable under pressure
- Short-setting characteristics and easy mixing saves valuable practice time







Superior hydrophilicity means superior performance!

In a recent study*, A-silicones and polyethers were investigated with respect to their wetting properties using video and high resolution contact angle measurements during the working time.

On the basis of initial hydrophilicity (quantified via the initial contact angle) and hydrophilic kinetics (spreading capacity of the water drop), both initial hydrophilic and high hydrophobic A-silicones were observed.

The time of initial contact and the speed with which the liquid spreads play a decisive role in improving the precision of an impression. Furthermore, it is important for the hydrophilicity to be retained throughout the entire working time of the material, not just at the beginning. The tested materials displayed a continuous decrease in hydrophilicity, while the hydrophilicity of the A-silicone Panasil[®] initial contact remained consistent throughout the entire working time.





* Comparative test carried out by Kettenbach contact angle measured using the measuring device

* Affinis Light™, Aquasil Ultra XLV™, and Impregum™ are not registered or trademarked by Kettenbach

DAS 100, from Kruss, 2006. Evaluated using the circle-fit method.

Comparison of spreading behavior begins 30 seconds after mixing starts and within the first 3 seconds after applying a water drop at 80% humidity. The test was performed on unpolymerized impression material.

Panasil[®] initial contact X-Light ▲ Contact angle 120° — Affinis™ light ■ Aquasil™ Ultra XLV 100° . Imprint™ 4 liaht **■** 80° 60° 40° 20° Panasil[®] initial contact X-Light Droplet age 0,5s 1s 2s 35

Measurements were carried out of the Kettenbach company on the premises on the basis of respective DIN/ISO standard measurement specifications.





Photos by Hubert Christian Roggendorf, D.D.S. – Department of Prosthodontics Center of Dento-Maxillo-Facial Medicine, Rheinische Friedrich-Wilhelms University of Bonn, Germany

Attains top values fastest

PANASIL®



Applications

4

		monophase	binetics Putty Fast	binetics Putty Soft	Putty	Putty Fast Set	Putty Soft	tray Fast Heavy	tray Soft Heavy
Dou	uble mix impression								
Fixa	ation impression								
Fun	nctional impression							-	
Two	o-step impression								
One	e-step impression								
	pression over fixed/ novable restorations	•	•		-		•		
	pression for forming ctional peripheries	-	-			•			•

highly recommended recommended

Technical Data

	Product name	Mixing ratio	Working time at 23° C/ 74° F ≤	Working time at 35° C/ 95° F ≤	Intraoral setting time at 35° C/ 95° F \geq	Total setting time* ≥	Hardness (approx.) Shore	Linear dimensional change (max.) %	Elastic Recovery test (approx.) %	Strain in compression (approx.) %
Mono-	Panasil [®] monophase	5:1	2 minutes	Not applicable	2 minutes	4 minutes	A 60	- 0.20	99.7	3.5
phase	Panasil [®] monophase	1:1	2 minutes	Not applicable	2 minutes	4 minutes	A 60	- 0.20	99.7	3.5
	Panasil [®] binetics Putty Fast	5:1	1 minute 30 seconds	Not applicable	2 minutes 30 seconds	4 minutes	A 63	- 0.20	99.5	2.5
	Panasil [®] binetics Putty Soft	5:1	2 minutes	Not applicable	3 minutes	5 minutes	A 56	- 0.20	99.5	3.5
Putty	Panasil [®] Putty	1:1	2 minutes	Not applicable	2 minutes	4 minutes	A 66	- 0.20	99.0	2.7
	Panasil [®] Putty Fast Set	1:1	1 minute 30 seconds	Not applicable	2 minutes	3 minutes 30 seconds	A 66	- 0.20	99.0	2.7
	Panasil [®] Putty Soft	1:1	2 minutes	Not applicable	2 minutes	4 minutes	A 60	- 0.20	99.0	2.7
	Panasil [®] tray Fast Heavy	5:1	1 minute 20 seconds	Not applicable	2 minutes	3 minutes 20 seconds	A 62	- 0.20	99.7	2.5
Tray/	Panasil [®] tray Fast Heavy	1:1	1 minute	Not applicable	2 minutes	3 minutes	A 62	- 0.20	99.5	3.0
Heavy	Panasil [®] tray Soft Heavy	5:1	2 minutes	Not applicable	2 minutes	4 minutes	A 55	- 0.20	99.7	3.0
	Panasil [®] tray Soft Heavy	1:1	2 minutes	Not applicable	2 minutes	4 minutes	A 55	- 0.20	99.5	3.0
	Panasil [®] tray Soft Heavy Fast	1:1	1 minute	Not applicable	2 minutes	3 minutes	A 55	- 0.20	99.5	3.0
	Panasil [®] initial contact Regular	1:1	1 minute 30 seconds	1 minute	2 minutes 30 seconds	4 minutes	A 46	- 0.20	99.7	3.0
	Panasil® initial contact Regular Fast	1:1	1 minute	30 seconds	2 minutes	3 minutes	A 46	- 0.20	99.7	3.0
Wash/	Panasil [®] initial contact Light	1:1	1 minute 30 seconds	1 minute	2 minutes 30 seconds	4 minutes	A 46	- 0.20	99.3	3.5
Light	Panasil [®] initial contact Light Fast	1:1	1 minute	30 seconds	2 minutes	3 minutes	A 46	- 0.20	99.3	3.5
	Panasil [®] initial contact X-Light	1:1	1 minute 30 seconds	1 minute	2 minutes 30 seconds	4 minutes	A 46	- 0.20	99.3	3.5
	Panasil [®] initial contact X-Light Fast	1:1	1 minute	30 seconds	2 minutes	3 minutes	A 46	- 0.20	99.3	3.5
Lab Putty	Panasil [®] lab Putty	1:1	1 minute 45 seconds	**	**	6 minutes	A 85	- 0.10	-	-
Fully	Panasil [®] lab Putty Soft	1:1	1 minute 45 seconds	**	**	6 minutes	A 73	- 0.10	-	-

* Total setting time (removal from the mouth) from the start of mixing. ** Laboratory product values not measured



Seven times maximum precision and comfort.

Fast-setting, syringeable bite registration material (A-silicone) with a high final hardness for exact recording of occlusal relationships.

Product properties

- Comfortable (Shore-A 62), hard (Shore-A 90) or extra-hard (Shore-D 43) material
- Quick, extra-quick or extra-long working time
- ✓ Short or extra-short intraoral setting time
- Elastomeric, stable storage, can be cut or trimmed, very thixotropic

Market leader in Germany¹⁾ over 20 years Futar[®] = over 20 years quality + over 20 years innovation

Futar® Clear Fast

Futar® Cut & Trim Fast

Futar® D Slov

5

User advantages

- ✓ Precise results, maximum convenience, top quality
- Product range provides the correct material for every indication
- No rebound when aligning the models

Benefits for the dentist

- Easy handling during application, removal and alignment
- Quick-setting characteristics save time in practice
- ✓ Wide range of products ensures the optimum material always
- Consistently top quality from the market leader = reliability of a precision bite registration



Technical Data

Product Information	Futar®	Futar [®] Fast	Futar [®] D	Futar® D Fast	Futar [®] D Slow	Futar [®] Cut & Trim Fast	Futar [®] Clear Fast
Working time at 23° C / 74° F \leq	30 seconds	15 seconds	30 seconds	15 seconds	1 minute 30 seconds	15 seconds	15 seconds
Intraoral setting time at 35° C /95° F ≥	1 minute 30 seconds	45 seconds	1 minute 30 seconds	45 seconds	3 minutes	45 seconds	45 seconds
Total setting time* \geq	2 minutes	1 minute	2 minutes	1 minute	4 minutes 30 seconds	1 minute	1 minute
Hardness (approx.) Shore	A 90	A 90	D 43	D 43	D 43	D 35	A 62

* Total setting time (removal from the mouth) from the start of mixing.

1) according GfK Germany

2) 5+ rating for Futar Fast only

Futar® F

utar® D Fas



The best of both worlds.

6

Heavy, medium and light body, elastomer impression material on a Vinylsiloxanether[®] basis for the one-step impression technique, combines the advantages of polyether and A-silicone.

Your real alternative to Polyether

Product properties

- The best of polyether
- ✓ The best of A-silicone
- ✓ Even shorter intraoral setting time

User advantages

- Pronounced flowability, optimum hydrophilicity
- Odorless and tasteless
- High final hardness immediately after setting ensures optimum fixation
- Extra-long intraoral working time, making it ideal for large restorations
- ✓ No gag reflex
- ✓ No distortion

Benefits for the dentist

- Excellent flowability which yields remarkable hydrophilicity
- ✓ Very short intraoral setting time saves valuable chair time
- Wide range of applications (pick-up, functional and crown/bridge impressions): you still only require one material!
- Easy removal from the mouth and model, ensures a stress-free impression-taking process

90 seconds quicker and more hydrophilic than Polyether!









Best hydrophilicity means best performance!

Smart hydrophilicity concept for perfect wetting



Thanks to its highly dynamic surface conditioning, Identium® reaches optimum values for hydrophilicity and offers exceptional wetting even in a moist environment with the lowest achievable contact angle (less than 10° after 1 second).

7

Initial contact angle



Comparison of spreading behavior begins 30 seconds after mixing starts within the first 3 seconds after applying a water drop at 80% humidity. The test was performed on unpolymerized impression material.

With a contact angle of less than 10° after one second, Identium® reaches maximum values for hydrophilicity. In contrast to other materials available in the market, Identium® succeeds in quickly developing high hydrophilicity and then maintaining it over the entire working time. The result: well-defined impressions, even in extreme situations.

The smaller the contact angle, the more hydrophilic is the behavior of the material.

Product name	Mixing ratio	Working time at 23° C/ 74° F ≤	Working time at 35° C/ 95° F ≤ 1)	Intraoral setting time at 35° C/ 95° F ≥	Total setting time* ≥	Hardness (approx.) Shore		Elastic- Recovery test (approx.) %	Strain in com- pression (approx.) %
ldentium® Medium	5:1 jumbo cartridge	2 minutes	80 seconds	2 minutes 30 seconds	4 minutes 30 seconds	A 60	- 0.20	99.0	2.3
Identium® Medium Fast	5:1 jumbo cartridge	1 minute 15 seconds	40 seconds	2 minutes 15 seconds	3 minutes 30 seconds	A 60	- 0.20	99.0	2.3
ldentium® Heavy	5:1 jumbo cartridge	2 minutes	80 seconds	2 minutes 30 seconds	4 minutes 30 seconds	A 60	- 0.20	99.0	2.8
Identium® Heavy Fast	5:1 jumbo cartridge	1 minute 15 seconds	40 seconds	2 minutes 15 seconds	3 minutes 30 seconds	A 60	- 0.20	99.0	2.8
ldentium® Light	1:1 cartridges	2 minutes	80 seconds	2 minutes 30 seconds	4 minutes 30 seconds	A 46	- 0.20	99.0	3.8
Identium® Light Fast	1:1 cartridges	1 minute 15 seconds	40 seconds	2 minutes 15 seconds	3 minutes 30 seconds	A 46	- 0.20	99.0	3.8

Technical Data

*Total setting time (removal time from mouth) from start of mix.

1) With a reduced working time the total setting time of 3 minutes 30 seconds for the fast set version and 4 minutes 30 seconds for the regular set version has always to be heeded.

8

VISALYS® CORE



Secure bonding with your adhesive.

Visalys[®] Core is the first core build-up material to feature the unique Active-Connect-Technology for bonding reliably to all popularly used brands of single and multi-stage adhesives. Therefore,

Visalys[®] Core even enables dual-curing core build-up materials to bond to universal adhesives – without an additional activator.

Using Visalys[®] Core for building up the core and cementing the root post creates a "seamless post and core" without weak points and consisting of only one material.

Top bond strength with all popularly used brands of adhesives!

Benefits for the dentist

- Active-Connect-Technology: Top bond strength with all popularly used brands of single and multi-stage adhesives – without an additional activator!
- \checkmark One material, two applications: Core build-up and cementation of root posts
- Excellent positional stability yet good flowability and low extrusion force
 Dual-curing: Saves time due to optional light-curing, reliable
- self-curing even in deep areas
- Highly stable due to very good compressive strength
- Minimal smear layer







The new standard for taking impressions of the opposing arch

Medium body, elastomer A-silicone with alginate-like consistency, specially developed for taking impressions for the opposing arch, temporary crowns/bridges and for the fabrication of splints.

Storable and can be poured several times

9

Product properties

- ✓ Optimum elastic recovery
- High dimensional stability
- ✓ Very thixotropic, high flowability and tear resistance
- ✓ Neutral or Strawberry scent

User advantages

- ✓ A-silicone with alginate-like properties
- Impression can be stored
- ✓ Maximum precision, easy removal from the mouth
- Easy handling with familiar components

Benefits for the dentist

- Bite registration, precision and opposing arch impression taking using the same material optimally coordinated
- Can be poured several times, even after a long period of storage
- Ideal impression material for the fabrication of temporary restorations and splints
- Standardized, hygienic procedures in line with quality management



For temporary crowns and bridges



Can be poured several times



VISALYS® TEMP



For strong temporary crowns and bridges.

Two-component material based on a multifunctional acrylic composite. Suitable for the fabrication of temporary crowns, partial crowns, inlays, onlays and veneers.

10

Product properties

- Very stable and fracture-resistant
- ✓ Setting with minimum inhibition layer
- ✓ Very smooth surface
- ✓ Free of Bisphenol A

Benefits for the dentist

- ✓ Long, problem-free time in situ, fewer repairs
- Easy handling, satisfied patients
- Easy to process, can be precisely trimmed with minimal dust formation
- No health exposure of patient and user; therefore meets the requirements of tomorrow.







Outstanding properties for very strong temporary restorations.

User advantages

- ✓ Also suitable for long-term temporary restorations
- Minimum smear layer, low heat build-up
- ✓ High inherent sheen, easy to polish
- No allergic potential due to softening agent

VISALYS® TEMP

KETTENB NTAL Simply intelligent

The result: exceptionally stable

High flexural strength combined with high impact strength reduces the risk of the temporary restoration fracturing in the mouth as a result of high or sudden masticatory loading -e.g. when biting unexpectedly on a cherry stone.

Impact strength	n [kJ/mm²]	
Visalys [®] Temp		12.5
Protemp [™] Plus		9.5
Luxatemp [®] Ultra		7.8
Structur [®] 3		8.5
Integrity®		6.8

Measurement at the University of Erlangen, Germany, according to DIN EN ISO 179, DIN 53453

Flexural strength [MPa]



Internal measurements Kettenbach, according to DIN EN ISO 10477 (Feed rate:1 mm/min)

lower the risk that, e.g. the margins of the temporary restoration may break off.

The higher the diametral tensile strength, the A high E-modulus means that the temporary restoration is resistant to deformation and reliably conserves the intraoral situation until the permanent restoration is fitted.

Diametral tensile strength [MPa]	
Visalys [®] Temp	67
Protemp [™] Plus	51
Luxatemp [®] Ultra	51
Structur [®] 3	47
Integrity®	49

Internal measurements Kettenbach, according to ADA Specification No. 27:76



Internal measurements Kettenbach, test set-up according to **DIN EN ISO 10477**

11

MUCOPREN® SOFT



Holds to the denture. And does what it promises.

Permanently soft reline material based on polyvinyl siloxane for direct and indirect application.

Product properties

- Smooth, hydrophobic silicone surface
- ✓ Silicone-based sealing
- ✓ High tear-resistance
- Patented adhesive

User advantages

- Resistant to contamination
- Good trimability
- Excellent bond strength to the denture
- For direct and indirect application

Benefits for the dentist

- ✓ Optimum requirements for hygienic denture conditions
- ✓ Long period in situ and comfort
- Modern material characteristics ensure easy, reliable handling
- ✓ Very efficient, value creation remains in the dental practice



Direct reline



Can be easily trimmed



Indirect reline



to seal



For denture relining without compromises!



Sugi[®] Butterfly – for Dental Tissue Isolation.

Uniquely prepared Versatile Cellulose Sponge Material made from a special formulation of pure cotton and regenerated cellulose, natural biocompatible components. Used for absorption of oral fluids and tissue isolation during dental procedures.

Product properties

- Highly absorbent (15 17 times)
- Made of pure cotton and regenerated cellulose
- ✓ Fast wicking
- ✓ No linting
- ✓ High tensile strength

User advantages

- ✓ Comfortable for the patient
- Remains soft and elastic
- Can be sterilized
- Biocompatible natural components

Benefits for the dentist

- ✓ Protects soft tissue during operation
- ✓ Keeps area dry, no suction needed
- ✓ Cost effective solution







The biocompatible

Sugi[®] Sponge has

been in medical use

for 50 years

AIRWAY METRICS -

System used to Obtain An Optimum Treatment Position for Oral Appliance Therapy.



- Quickly and Easily Determine an Initial Treatment Position
- Multiple Horizontal and Vertical Positions
- Designed Specifically for Dental Sleep Medicine

For more information contact Kettenbach direct 877-532-2123







www.kettenbachusa.com

Product	REF	piece pack
Snore Screener/ Mandibular Positioning Simulator System Blue	REF 65100	16 piece pack
Snore Screener/ Mandibular Positioning Simulator System Amber	REF 65101	16 piece pack
Snore Screener/ Mandibular Positioning Simulator System Black	REF 65102	16 piece pack
Bite Fork Pack & Handles	REF 65103	24 piece pack (23 regular, 1 wide, 4 handles)
Scanning Bite Fork Pack & Handles	REF 65104	15 piece pack (Includes 4 handles)
Vertical Titration Key Blue	REF 65105	9 piece pack
Vertical Titration Key Amber	REF 65106	9 piece pack
Vertical Titration Key Black	REF 65107	9 piece pack

Product	1 x 5 ml double syringe, 10 brown mixing tips, 5 intraoral tips,	2 x 5 ml double syringe, 20 brown mixing tips, 10 intraoral tips,	Normal pack: 25 ml, 1:1, 1 x 25 ml cartridge, 20 yellow mixing tips, 20 intraoral tips
Visalys [®] Core, white	-	REF 13860	REF 13870
Visalys® Core, dentin	-	REF 13861	REF 13871

ASSORTMENT

Visalys® Temp, shade BL



Intro pack, Regular set	Refill pack, Regular set	Intro pack, Fast set	Refill pack, Fast set		
380 ml + accessories	2 x 380 ml	380 ml + accessories	2 x 380 ml		
REF 14712	REF 14713	- -	-		
REF 14714	REF 14715	-	-		
REF 14716	REF 14717	REF 14718	REF 14719		
REF 14724	REF 14725	REF 14726	REF 14727		
REF 14708	REF 14709	-	-		
-	-	REF 14700	REF 14701		
REF 14702	REF 14703	-	-		
-	-	REF 14704	REF 14705		
REF 14706	REF 14707	-	-		
Normal pack: 900 ml		Economy pack			
REF	11101		-		
REF	11141		-		
REF	11121		-		
	-	REF 1	1153 (2 x 5 kg)		
	-	REF	1155 (2 x 5 kg)		
A	Normal pack		Bulk pack,		
	Regular set (2 x 50 ml)		r set (24 x 50 ml)		
-	REF 13501		-		
-	REF 13561		REF 13563		
-	REF 13551				
-	REF 13541				
REF 13846 (6 x 50 ml)	-				
REF 13826 (6 x 50 ml)	-		REF 13827		
-	REF 13701/REF 13711		-		
-	REF 13401/REF13461	REF 1	3402/REF 13463		
-	REF 13411/REF13471	REF 1	3412/REF 13473		
-	REF 13431/REF13481		-		
-	REF 11912		-		
-	REF 11926		REF 11927		
-	REF 11932		-		
-	REF 11961		REF 11967		
-	REF 11951		-		
-	REF 11975		-		
-	REF 11928		-		
-	REF 15687	-			
REF 28105 – 50 ml M (Ø 6.5 mm), 20 blue	lucopren Soft, 50 ml Mucopren sil e mixing tips (Ø 3.2 mm), 1 brush	licone sealant, 10 ml Mucoprer holder, 20 single-use brushes,	adhesive, 7 green mixing tips 1 steel cutter, accessories		
	RE	F 13780			
	RE	F 13781			
	nL	1 10/02			
		F 13790			
	380 ml + accessories REF 14712 REF 14714 REF 14714 REF 14716 REF 14708 - REF 14706 Normal pack: 900 ml REF REF <	380 ml + accessories 2 x 380 ml REF 14712 REF 14713 REF 14712 REF 14713 REF 14714 REF 14715 REF 14716 REF 14717 REF 14708 REF 14709 - - REF 14708 REF 14709 - - REF 14702 REF 14703 - - REF 14706 REF 14707 Normal pack: 900 ml - Normal pack: 900 ml REF 11101 REF 11121 - - - REF 1300 REF 13501 REF 13551 - - - - - - - - - - - - - REF 13826 (6 x 50 ml) - - REF 13846 (6 x 50 ml) - REF 13841/REF 13711 REF 13826 (6 x 50 ml) - - REF 13401/REF 13711 - REF 13841	380 ml + accessories 2 x 380 ml 380 ml + accessories REF 14712 REF 14713 - REF 14712 REF 14713 - REF 14714 REF 14715 - REF 14716 REF 14717 REF 14718 REF 14708 REF 14709 - REF 14708 REF 14709 - REF 14700 REF 14703 - REF 14702 REF 14703 - REF 14706 REF 14703 - REF 14706 REF 14707 - Normal pack: 900 ml Economy pack Economy pack REF 11121 - REF 11121 REF 11121 - REF 11121 REF 11121 - REF 11121 REF 13551 - REF 11121 REF 13551 - REF 11121 REF 13846 (6 x 50 ml) - REF 13551 REF 13846 (6 x 50 ml) - - REF 13846 (6 x 50 ml) - - REF 13846 (6 x 50 ml) - - REF		

REF 13788

Kettenbach LP 16052 Beach Blvd. Ste. 218 Huntington Beach, CA 92647, USA Phone: 877-532-2123 Fax: 888-892-9820 Custservice@kettenbachusa.com www.kettenbachusa.com



022763_3019 USA