EN Esthetic restorations in severely darkened anterior teeth using Whitepost and Brava Block.

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QUEJA PRINCIPAL

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INITIAL CLINICAL / X-RAY EVALUATION

Tooth 11 presented a severe degree of darkening, poor composite restoration but with adequate endodontic treatment. The gingival margins of the central incisors were misaligned and the clinical crowns of all anterior teeth were short.

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TREATMENT PERFORMED

A crown increase was planned from element 13 to 23 along with the regularization of the gingival margins and cementation of an intraradicular fiberglass retainer. In addition, a total crown was made of glass-ceramic composite for CAD/CAM (Brava Block - FGM) for tooth 11.

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INTRODUCTION

Darkened teeth are always a challenge in the clinical routine and, when deciding the type of treatment, one should always select the most conservative one. In this context, whitening becomes an interesting alternative. However, due to the great relapse of this type of treatment, the patient usually ends up choosing a more definite solution such as crowns and veneers.

Fig. 1 - Initial smile.
Fig. 1 - Sonrisa inicial.

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CASE REPORT

A female patient attended the subject of Integrated Clinic of the Fluminense Federal University (Rio de Janeiro, Brazil), with the request for the realization of an anterior crown with the appropriate format for future orthodontic finalization.

Endodontically treated teeth are susceptible to biomechanical failures, especially when there is a significant loss of dental structure¹. In such cases, the use of intraradicular retainers is necessary, mainly because the remnant structure is insufficient to adequately support and retain the final restoration^{2,3}. Several techniques are available for the restoration of teeth with endodontic treatment and, within this context, choosing an inadequate method when restoring these teeth may compromise their longevity.⁴

The tooth to be restored presented a good amount □

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► of coronal structure and conservative endodontic treatment, therefore the direct technique was indicated (see table 1).

The conduit was cleaned and dried with paper cones for the adhesive procedure and light-curing was performed for 40 seconds. With the aid of an insertion tip, Allcem Core (FGM) dual cement was applied in the conduit from the apical region to the coronary chamber. The post Whitepost DC-E 2 (FGM) was inserted and the set was light cured for 40 seconds. The excess of the post was cut and the occlusal adjustments were performed. The patient was submitted to periodontal surgery by Professor Ronaldo Barcellos at the Fluminense Federal University (Rio de Janeiro – Brazil). After 90 days, total crown preparation and molding were performed and sent to the laboratory. A total crown was made of Brava Block (FGM), glass-ceramic composite for CAD/CAM, with external characterization. After tests and adjustments, the crown was cemented.

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Fig. 4 - Removal of gutta-percha with Peeso drill number 3, maintaining 5mm in the apical region for safety of the endodontic treatment.

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Fig. 5a and 5b - Conduit adequacy with Whitepost DC-E 2 (FGM) drill for the perfect juxtaposition of the fiberglass post.

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Clinical Tip:

Scaling of drills provides a gradual wear of structure and facilitates the preparation of the conduit.

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 German technology.
Tecnología alemana.

Fig. 6a and 6b - The DC-E 2 post was placed in position. The excellent adaptation and complete filling of the conduit can be observed.

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Whitepost


Adapts perfectly: to the conduit, to your technique and to you.

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Postes de fibra de vidrio.

Prosil

Increases adhesion to different resinous substrates.

Silanizing agent for ceramics, ceromers, laboratory composites and fiberglass posts.

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Fig. 7 - The post was cleaned with 70% alcohol and silanized (Prosil - FGM).

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Fig. 8a and 8b - The conduit was cleaned and dried with paper cones for the adhesive procedure.

Fig. 9 - Ambar Universal APS (FGM) adhesive was applied actively inside the conduit in the self-etching mode, the excesses removed with the aid of a Cavibrush (FGM) and then light-cured for 40 seconds.

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Fig. 10 - With the aid of a tip, Allcem Core (FGM) dual cement was applied in the conduit from the apical region to the coronary chamber.

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Allcem Core

Triple proven efficiency.

Dual resin cement for core build-up and adhesive cementation of intraradicular posts and prosthetic crowns.

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Fig. 11 - DC-E 2 post was inserted into the conduit filled with cement.

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Fig. 12a and 12b - The set was light-cured for 40 seconds. The excess of the post was cut and the occlusal adjustments were performed.

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Fig. 14 - After preparation and the provisional piece with the new format were carried out, a small adjustment in the gingival margin of tooth 21 was performed with electric scalpel.

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FGM MATERIALS USED / MATERIALES FGM UTILIZADOS:

REFERENCES / REFERENCIAS

1. Ferro, M.C.L. et al. Fracture strength of weakened anterior teeth associated to different reconstructive techniques. *Brazilian Dental Journal* 2016; 27(5):556-561.
2. Gulerenler, K.A. et al. Long term clinical outcomes of endodontically treated teeth restored with or without fiber post: retained single unit restorations. *Journal of endodontics*, 2017;43(2).
3. Clavijo,V.G.R. et al. Fracture strength of flared bovine roots restored with different intraradicular posts. *Journal of applied oral Science*. 2009;17(6):574-8.
4. Santos-Filho, P.C.F. et al. Influence of ferule, post system, and lenght on biomechanical behavior on endodontically treated anterior teeth. *JoE* 2014; 40(1).